

# VAV INSTRUMENTS (SAFE AND SECURE)

## AUTOPSY TABLE (POSTPARTUM)



The complete autopsy workstation is made of high grade Stainless Steel. The complete table is aesthetically designed by using state-of-the-art engineering, heavy gauge stainless steel and conveniences. The stainless steel table is mounted on stainless steel tubular stand. The autopsy table is scientifically designed for maximum efficiency.

Approximate Overall Dimensions:  
85"Length X 36"Width X 36"Height

(Top-to-bottom Heavy Gauge Stainless Steel)Eliminating iron and paint in the construction means a more durable, longer-lasting finish.

The top-to-bottom stainless steel construction will withstand the harshest conditions.

Sink is integrally constructed of 18-gauge stainless steel with rounded corners for easy cleaning.

Work surfaces are sloped to facilitate drainage into the extra large integral sink.

All switches are encased in a way to assure easy clean up.

We install the highest quality plumbing fixtures for innovative design, quality craftsmanship and cutting-edge safety features.

**Power Requirement:** - The instrument is designed to operate on 220V AC Supply