

BLOOD BAG TUBE SEALER



BLOOD BAG TUBE SEALER IS A COMPACT INSTRUMENT TO SEAL THE BLOOD BAG PILOT TUBE BY RADIO FREQUENCY SEALING SYSTEM, NO HAEMOLYSIS AT BLOOD IN THE TUBE BECAUSE RADIO FREQUENCY SEALING SEALS THE TUBE WITHOUT DAMAGING THE BLOOD CELLS INSIDE, AND EASY SEPARATION OF TUBE THE SEALING IS DONE, SEALING CAN BE ACHIEVED IN LESS THAN 1.5 SEC. PER TUBE BY SIMPLE OPERATION.

Name of the equipments	Blood bag Tube Sealer
Oscillation Frequency	40.68 MHz
Radio Frequency output	40.68 MHz
Sealing time	Approx 1.5 sec.
Input Volt	230 V AC
Size	180mmX215mmX245mm(HXD _X W)