

PLATELET AGITATOR WITH INCUBATOR.



PROVIDES AN IDEAL STORING ENVIRONMENT FOR THE PLATELETS. THIS UNIT DESIGNED TO MAINTAIN A STABLE AND UNIFORM TEMPERATURE OF 22.0°C, THE USER-FRIENDLY DESIGN INCLUDES AN ACRYLIC INNER DOOR, AND SEE THROUGH OUTER DOOR, SIMPLE CALIBRATION IS POSSIBLE VIA THE CONTROL PANEL. THIS, COUPLED WITH A VISUAL AND AUDIBLE ALARM SYSTEM, ENSURES THAT YOUR SAMPLES ARE BOTH SAFE AND SECURE, BUT ALSO EASY TO ACCESS WHEN YOU NEED THEM

Other features of this unit include :

Micro Controlled Temperature controller with Audio Visual Alarm for high & Low Temperature.

Display of Set and Process value

Extra-strength non CFC refrigeration System

The unit is mounted on castor wheel for easy mobility.

CFC Free PUF (Poly Urethane Foam) Insulation to ensure temperature stability and save power consumption

Inner Dim(WxDxH)	Operating Temperature
18x18x16inch .	22.0 °C (48 PLATELET BAGS)