



13485:2016



## About Us

Established in the year 2006 at Maharashtra, India, we “VAV Instruments” are Sole Proprietorship (Individual) based firm, involved as the Manufacturer and Exporters of Blood Bank Instruments manufacturer in Mumbai, Laboratory Freezers, Mortuary Chamber manufacturer in Mumbai India, Ice Lined Refrigerator Manufacturer in Mumbai India, Blood Donor Couch, Dead Body Freezer Manufacturer in Mumbai, Laboratory Incubators, Blood Collection Monitor, OT Table, Blood Bank Scale, Circulating Water Bath etc. All our products are getting widely acclaimed among the large clientele for their exclusive designs, superior quality, and reliability. Apart from this, our ability to maintain timelines as well as quality in the assortment, providing cost-effective solutions and assurance to make timely shipment of the orders placed by customers have assisted us positioning our name in the list of top-notch companies of the industry.



## ICE LINED REFRIGERATORS

Ice lined refrigerators are designed for storing life saving vaccines. our sturdy and reliable Ice Lined refrigerators with unique design innovations are specially designed to store vaccine products at a controlled, stable temperature and provides a means of storing vaccines and medical supplies. These Ice Lined Refrigerators with a Temperature range of 2 Deg. C to 8 Deg. C come with a highly efficient backup system of 24 – 72 hours. our Ice Lined Refrigerators are designed to maintain the COLD CHAIN which is a must feature of Ice Lined refrigerators for routine & emergency vaccine handling.

Strong and reliable refrigeration is the main feature of ice lined refrigerator. The unit is fitted with a digital electronic temperature controller cum indicator for easy readability along with high & low alarm system. The heart of the system is a highly efficient compressor and condenser with state of the art accurate Temperature Controls constrained to work with low energy consumption.

### FEATURE SELECTIONS:

1. High and low temperature alarms.
2. Temperature data logger.
3. Temperature Chart Recorder.
4. Computer Interface & software.
5. Data store memory & print function.
6. Alarm system with rechargeable battery backup.
7. Alarms for open door, sensor fault, condenser fault etc.
8. Voltage safety system.
9. Dual door system with inner Glass door.
10. Vacuum insulated Glass door system.

MODELS	VILR60	VILR150	VILR 200	VILR250	VILR300	VILR350	VILR400
TEMPERATURE	2°C TO 8°C	2°C TO 8°C	2°C TO 8°C	2°C TO 8°C	2°C TO 8°C	2°C TO 8°C	2°C TO 8°C
ACCURACY (- / +)	5%	5%	5%	5%	5%	5%	5%
INNER SIZE CM (W/D/H) APP.	44/47/42	70/50/50	80/60/45	90/60/50	90/70/50	115/65/45	117/70/50
BACKUP TIME(IN HRS)	24 TO 72	24 TO 72	24 TO 72	24 TO 72	24 TO 72	24 TO 72	24 TO 72
CAPACITY IN LTR. (APP.)	86 liter	175 liter	216 liter	270 liter	315 liter	362 liter	409 liter

